



Product Data Sheet

AMBERLITE™ FPA90 CI Ion Exchange Resin

Food- and Biopharmaceutical-grade, Styrenic, Macroporous, Strong Base Anion Exchange Resin for Cane Sugar Decolorization and Biopharmaceutical Processing

Description

AMBERLITE™ FPA90 CI Ion Exchange Resin has been specially designed for the decolorization of liquid sugar syrups. Ion exchange based decolorization technology has been used more effectively and economically than carbon or bone char based technologies. Sugar refiners and soft drink bottlers around the world have installed AMBERLITE™ FPA90 CI to successfully decolorize sucrose solutions.

AMBERLITE™ FPA90 CI is a styrenic, macroporous, Type I strong base anion resin.

Cane Sugar Decolorization

AMBERLITE™ FPA90 CI Resin is an excellent choice for cane sugar decolorization, offering advantages such as:

- A leading resin used in cane decolorization, with many years of successful service
- Exceptional physical stability, excellent resistance to osmotic shock, and very good organic fouling resistance

Styrenic AMBERLITE™ FPA90 CI can be used as a single component or following acrylic AMBERLITE™ FPA98 CI Ion Exchange Resin for highly-colored feed solutions.

Biopharmaceutical Processing

AMBERLITE™ FPA90 CI Resin is an excellent resin of choice for decolorization of high molecular weight organic color bodies in many bioprocessing applications such as natural product extraction and recovery of antibiotics from fermentation broth. It is commonly used in aminoglycoside purification bioprocess as well as in macrolide antibiotics processes like erythromycin, and tylosin, the latter being primarily used in animal health. Derivatives of erythromycin, including clarithromycin and azithromycin, are particularly useful in treating respiratory infections.

A combination of AMBERCHROM™ CG Chromatography Resins and/or AMBERLITE™ FPC3500 Ion Exchange Resin with AMBERLITE™ FPA90 CI allows a higher antibiotic purity level due to the decolorization capabilities of the latter either in the pre- or post-purification step (vancomycin broth decolorization).

Applications

- Cane sugar decolorization
- Antibiotic decolorization

Typical Properties

| Physical Properties | |
|----------------------------|--|
| Copolymer | Styrene-divinylbenzene |
| Matrix | Macroporous |
| Type | Strong base anion, Type I |
| Functional Group | Trimethylammonium |
| Physical Form | Off-white, opaque, spherical beads |
| Chemical Properties | |
| Ionic Form as Shipped | Cl ⁻ |
| Total Exchange Capacity | ≥ 1.00 eq/L |
| Water Retention Capacity | 58 – 64% |
| Particle Size [§] | |
| Particle Diameter | 650 – 820 μm |
| < 300 μm | ≤ 0.5% |
| > 1180 μm | ≤ 3.0% |
| Stability | |
| Swelling | Cl ⁻ → OH ⁻ : ~25% |
| Density | |
| Particle Density | 1.050 – 1.080 g/mL |
| Shipping Weight | 700 g/L |

[§] For additional particle size information, please refer to the [Particle Size Distribution Cross Reference Chart](#) (Form No. 177-01775).

Suggested Operating Conditions

| | |
|--|---|
| Maximum Operating Temperature (Cl ⁻ form) | 80°C (176°F) |
| Bed Depth, min. | 1000 mm (3.3 ft) |
| Flowrates | |
| Service | 2 – 4 BV*/h (or up to 8 BV/h) |
| Sweeten-off | Service flowrate for 1.5 – 2 BV |
| Backwash | See Figure 1 |
| Regeneration | 2 BV/h |
| Slow Rinse | 2 BV/h |
| Sweeten-on | Service flowrate for 1 BV |
| Contact Time | |
| Regeneration | ≥ 45 – 60 minutes |
| Displacement Rinse | ≥ 60 minutes |
| Total Rinse Requirement | 5 BV |
| Regenerant | |
| Concentration | 10% NaCl 0.2% NaOH |
| Level, 100% basis | |
| Co-current | 180 – 200 kg/m ³ (11.3 – 12.5 lb/ft ³) |
| Counter-current | 150 kg/m ³ (9.4 lb/ft ³) |
| Temperature | 25 – 70°C (77 – 158°F) |

* 1 BV (Bed Volume) = 1 m³ solution per m³ resin or 7.5 gal solution per ft³ resin

Refer to the brochure [Ion Exchange Resins for Cane Sugar Decolorization](#) (Form No. 177-03556) for additional information.

Hydraulic Characteristics

Bed expansion of AMBERLITE™ FPA90 CI Ion Exchange Resin as a function of backwash flowrate and temperature is shown in Figure 1.

Pressure drop data for AMBERLITE™ FPA90 CI as a function of service flowrate and viscosity is shown in Figure 2. These pressure drop expectations are valid at the start of the service run with clean feed and a well-classified bed.

Figure 1: Backwash Expansion
Temperature = 5 – 40°C (41 – 104°F)





Figure 2: Pressure Drop
Viscosity = 1 – 10 cP



Product Stewardship

DuPont has a fundamental concern for all who make, distribute, and use its products, and for the environment in which we live. This concern is the basis for our product stewardship philosophy by which we assess the safety, health, and environmental information on our products and then take appropriate steps to protect employee and public health and our environment. The success of our product stewardship program rests with each and every individual involved with DuPont products—from the initial concept and research, to manufacture, use, sale, disposal, and recycle of each product.

Customer Notice

DuPont strongly encourages its customers to review both their manufacturing processes and their applications of DuPont products from the standpoint of human health and environmental quality to ensure that DuPont products are not used in ways for which they are not intended or tested. DuPont personnel are available to answer your questions and to provide reasonable technical support. DuPont product literature, including safety data sheets, should be consulted prior to use of DuPont products. Current safety data sheets are available from DuPont.

Please be aware of the following:

- **WARNING:** Oxidizing agents such as nitric acid attack organic ion exchange resins under certain conditions. This could lead to anything from slight resin degradation to a violent exothermic reaction (explosion). Before using strong oxidizing agents, consult sources knowledgeable in handling such materials.

Have a question? Contact us at:

www.dupont.com/water/contact-us

All information set forth herein is for informational purposes only. This information is general information and may differ from that based on actual conditions. Customer is responsible for determining whether products and the information in this document are appropriate for Customer's use and for ensuring that Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. The product shown in this literature may not be available for sale and/or available in all geographies where DuPont is represented. The claims made may not have been approved for use in all countries. Please note that physical properties may vary depending on certain conditions and while operating conditions stated in this document are intended to lengthen product lifespan and/or improve product performance, it will ultimately depend on actual circumstances and is in no event a guarantee of achieving any specific results. DuPont assumes no obligation or liability for the information in this document. References to "DuPont" or the "Company" mean the DuPont legal entity selling the products to Customer unless otherwise expressly noted. NO WARRANTIES ARE GIVEN; ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED. No freedom from infringement of any patent or trademark owned by DuPont or others is to be inferred.

DuPont™, the DuPont Oval Logo, and all trademarks and service marks denoted with ™, ℠ or ® are owned by affiliates of DuPont de Nemours Inc. unless otherwise noted. © 2019 DuPont.

